

EASY WALK AFO-10E eco



Easy Walk® AFO-10E is similar to AFO-10A version, but produced by combining carbon fibre with composite materials. Therefore it is the most cost-efficient from the Easy Walk® line. However, comparison tests have confirmed that the product keeps the functionality and comfort level close to pure carbon fibre version, with the exception of slight difference in weight. Comparing with pure carbon orthoses, due to lower resilience of composite materials, the *spring* effect is slightly lighter, what may serve as an advantage in certain cases.

Indications:

- drop foot secondary to neurological pathologies;
- foot and calf extensors after operation;
- ankle arthritis.

Contraindications:

- calf varicose;
- inflammatory lesions of calf skin;
- oedema;
- severe ankle and foot deformities;
- mental diseases.

Article Number	Side	Size	Shoe size	Height (cm)
AFO10E=LXS	Left (L)	XS	34-36	33
AFO10E=LS	Left (L)	S	36-39	35
AFO10E=LM	Left (L)	M	39-42	36
AFO10E=LL	Left (L)	L	42-45	37
AFO10E=LXL	Left (L)	XL	45-47	39
AFO10E=RXS	Right (R)	XS	34-36	33
AFO10E=RS	Right (R)	S	36-39	35
AFO10E=RM	Right (R)	M	39-42	36
AFO10E=RL	Right (R)	L	42-45	37
AFO10E=RXL	Right (R)	XL	45-47	39